


Example Comparison Report

Clean no-finding case

Case: 302b0d070d5150b91bafc935e94a847b

Review status: supported Confidence: high

Image Comparison



Model Comparison

Reference No finding	Swin Predicted: No finding Correct: No finding Missed: none	YOLO Predicted: No finding Correct: No finding Missed: none
--------------------------------	---	---

Findings

- Both the Swin classifier and YOLO detector demonstrate strong agreement in identifying no acute abnormalities on this chest radiograph.
- The Swin classifier shows a very high probability (99.
- 87%) for 'No finding' with all other pathological findings having very low probabilities (all <0.
- 05%).

Impression

- No acute cardiopulmonary abnormality identified on this chest radiograph.
- Both AI models demonstrate strong concordance in this assessment.

Agentic Workflow Summary

The agentic report favors No finding with a high confidence band and a supported review status. Primary AI-supported findings are none, while the remaining uncertain findings are none. Relative to radiologist labels, Swin captured No finding and YOLO localized No finding. Overall pattern emphasis is No acute abnormality.

Supporting Evidence

Top Swin Probabilities

- - No finding: 0.999
- - Nodule/Mass: 0.000
- - Pulmonary fibrosis: 0.000
- - Lung Opacity: 0.000
- - Other lesion: 0.000
- - Calcification: 0.000

Top YOLO Detections

- - None

Takeaways

- Swin correct labels: 1 / 1
- YOLO correct labels: 1 / 1
- Swin key miss: none
- YOLO key miss: none

Confidence Rules Used

- Swin positive label cutoff: 0.50
- YOLO detection confidence cutoff: 0.25
- Agentic AI confidence band is qualitative and selected by the reviewer.

This is an AI-generated decision-support summary that requires human radiologist review and interpretation for clinical use. This analysis should not be used as a substitute for professional medical judgment or definitive diagnosis.